

How to use your Easibreathe

- Shake the device well
- Hold the inhaler upright and open the cap
- Take care not to cover the air vent on the top of the device
- Take a breath in then breathe out fully.
- As you place the mouthpiece in your mouth, take a long, strong, deep breath in for as long as you can. The device will release a dose of medication in response to your breath - It is important that you continue to breath in when this happens.
- Remove from your mouth, holding your breath for up to 10 seconds.
- Close the cap – the device wont work if you don't
- Wait for 30 seconds then repeat until the prescribed number of doses have been taken

